Phosphate Buffered Saline with Tween® 20 (PBST-20X)

Supplied as 20X stock **✓** 1000 ml



Orders 877-616-CELL (2355)

orders@cellsignal.com

Support 877-678-TECH (8324)

info@cellsignal.com

Web www.cellsignal.com

rev. 11/07/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Phosphate Buffered Saline (PBS) solution with the detergent Tween® 20 for use as a wash buffer and diluent for ELISA. 1 X PBST contains 3.2 mM Na, HPO, , 0.5 mM KH₂PO₄, 1.3 mM KCl, 135 mM NaCl, 0.05% Tween® 20, pH 7.4.

Directions for Use: Dilute Phosphate Buffered Saline with Tween® 20 (PBST - 20X) to a 1X solution using ddH₂O. This product supplies enough 20X material to make 20 liters of

Storage: Supplied as a 20X solution. Store at room

Please visit www.cellsignal.com for a complete listing of recommended companion products.

© 2013 Cell Signaling Technology, Inc. Cell Signaling Technology® is a trademark of Cell Signaling Technology, Inc.